


DNA Ploidy Analysis



Ariol[®]




Automated scanning, segmentation,
and quantification of the DNA content
of tumor cells stained by Feulgen
reaction

Powerful DNA Ploidy Analysis

Beyond accurately measuring DNA Ploidy, Ariol provides additional controls to make analysis easier:

- Easy DNA quantification at 40X magnification
- Fast and automatic segmentation of nuclear boundaries
- Easy selection of tumor and control areas
- Split and join cells manually to enhance accuracy
- Sort DNA data by cell size and shape parameters
- Set thresholds to gate cells in and out by size and shape

DNA Ploidy Analysis



A real-time histogram shows the IOD of control cells overlaid with test cells. Interactive graphing tools can be used for viewing and selecting cells.



Fast and easy to use. The essential partner in DNA ploidy analysis.



The robotic arm of the SL-50 slide loader provides true automation

Example Workflow

- System retrieves case information via barcode, loads slide onto stage, pre-scans at 1.25x to locate tissue, and then auto-scans at 5x
- User selects the region of interest to analyze
- System auto-scans selected area at 40x
- User identifies areas for control tissue and tumor tissue
- System auto-segments nuclear boundaries within the selection areas, calculates the respective DNA indexes, and plots the case histogram
- User reviews the data to accept the nuclear segmentation as is or fine tunes it by splitting, joining, or redefining selective nuclear boundaries
- System automatically recalculates the DNA index and adjusts the histogram
System generates report

Time saving steps when system runs unattended User interactive steps

North America

Applied Imaging Corp.
120 Baytech Drive
San Jose, CA 95134-2302
USA

Toll-free : +1 800 634 3622

Telephone: +1 408 719 6400

Fax: +1 408 719 6401

International

Applied Imaging International Ltd
BioScience Centre
Times Square
Newcastle Upon Tyne NE1 4EP
UK

Telephone: +44 (0) 191 202 3100

Fax: +44 (0) 191 202 3101

