

The scanner, foundation of a good digital pathology solution

Ferenc Szipőcs
3DHISTECH - Sysmex



A good digital pathology solution

Scanner

- Fast
- Excellent image quality
- Compact, BF and FL in one package

Server

- Instant availability on LAN, WAN
- Management and Organization of slides
- Supports pathology workflow
- Integrates with HIS/LIS system
- Pathology workstation
- Image analysis
 - Quantification
 - Report




How soon can you see a new slide?

- Pannoramic 250 Super Flash
 - View a new slide in 100 seconds on average @
 40x
 - Winner of both 20x and 40X high throughput categories

Scan only once in 40x



Real time Scanning videos



Pannoramic Flash

Pannoramic Super Flash

Important to get high quality images

- Winner in the 40X category = 0.25 µm/pixel
 - Pannoramic DESK

- Zeiss Plan-Apochromat objectives
 - High Numerical aperture of 0.8 and 0.95
- Optics are the same in all of our products




Area sensor vs. TDI sensor

Area scanners

Camera sensor (1600x1200 pixel)

Line scanners


(4096 pixel) (≥64 x 4096 pixel)






Area camera advantages


- Very good co-localization
 - Important for FISH



- Z-Stack, BF and FL
 - More than 29 levels
 - 0.2 micron min step size
- Extended focus









Dedicated cameras for BF and FL

- Allows for excellent image quality
- Allows for high dynamic range
- Automated switch in less than 2 seconds







Fluorescence

- **✓9 Arbitrary FL Filters**
- **✓ Lumencor light source**
 - Solid-State, 6 channels
 - brighter than other lamps
- ✓ PCO Edge sCMOS camera
 - 16-bit, 2.5 megapixel, 40 fps
 - low noise, controls Lumencor
- **✓FISH**









16-bit Fluorescence imaging


Competitors have only max 10-bit

 Capture faint signals without washing out strong ones – high dynamic range

- Allows for high quality image processing without loosing relevant information
- JPEG 2000



Improved Image compensation in FL





Dark Field Preview

- Makes FL stained samples detectable
- Reduces false signals from fingerprints, dust, etc.
- Speeds up FL acquisition
 - We know where the sample is!
- Patented









Compact Scanners for every need



Pannoramic DESK
1 slide
BF only



Pannoramic MIDI 12 Slides BF / FL



Pannoramic 250 FLASH 250 Slides BF / FL



- Development since 1996
- 4th generation scanners
- Market leader in Europe
- Second in worldwide sales
 - Over 500 units sold
- Best Components, Expertise
- Collaboration with pathologists
- Workflow support
 - HIS/LIS integration



Thank you for your attention!

Questions?

Please visit the Sysmex booth for live scans and demos!

http://www.3dhistech.com

http://pathonet.org

http://slides.3dhistech.com/casecenter

http://www.sysmex-lifescience.com/

http://scanner-contest.charite.de/en/